


Bio-logic Navigator Pro AEP – designed by audiologists for audiologists.

One small package meeting your needs...



Multiple protocols available or create your own

Dual, scrolling displays for unlimited waveforms



Advanced data collection and analysis features

- Automatic electrode switching – eliminating the need to manually switch electrodes during testing
- Residual noise calculation and FsP – optimizing response detection
- Waveforms saved automatically – no lost data
- Flexible data analysis – mark and compare waveforms during collection
- Collection protocols – select from several built in protocols with commonly used parameters or easily create your own
- Quantifiable reproducibility – mathematically correlate waveforms to ensure reliability
- GraphMaster – create charts to compare latencies to published age referenced normative data




Flexible report software

- Save reports as pdf file – allowing for simple transfer to the patient's electronic medical record (EMR)
- Customizable report templates – creating professional reports to meet the needs of your practice
- Comprehensive patient database – allowing for creation of optional folders to organize patient records

Full AEP functionality includes:

- Click evoked ABR
- Tone burst ABR
- Electocochleography (ECoChG)*
- Middle and late latency response (MLR, ALR)
- P300
- Electrical ABR (EABR)*
- Vestibular Evoked Myogenic Potential (VEMP)**

*required cables available at additional cost
 **not available in the U.S.



The Bio-logic Navigator Pro family of products

AEP

Standard auditory evoked potentials including ABR, ECochG, MLR, ALR, P300 and EABR

Scout®

Full range of OAE test capability including DPOAE and TEOAE

ABaer®

Newborn hearing screening system including ABR, DPOAE and TEOAE

MASTER II

Multiple Steady-State Evoked Response (ASSR) technology

ENoG

External SEP stimulator option supports electroneurography

VEMP*

EMG rectification for Vestibular Evoked Myogenic Potential test

*not available in the U.S.



Scout



ABaer



MASTER II



ENoG



VEMP

Distributed & Marketed by



Bio-logic® Navigator® Pro

auditory evoked potentials



The **Bio-logic Navigator Pro AEP** is flexible, portable and adaptable, providing results you can trust.

Designed by audiologists to grow with you and your patients' needs

The Bio-logic Navigator Pro offers a full range of AEP functionality.

- Portable and light weight
- Multiple waveform display configurations
- Active data analysis during collection
- Continuous EEG display
- Customizable protocols and reports
- Windows® 7 compatible



GenWorks HEALTH PRIVATE LIMITED
Unit 522-524, 5th Floor, Gamma Block,
Sigma Tech Park, Bangalore - 560 066, India
Contact : Chandrakant Bangre
Mobile : +91 83191 54706